

Product datasheet for **TA309099**

Gja1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to amino acid residues from the C-terminal region conjugated to KLH
Formulation:	100 µl in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml BSA and 50% glycerol.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	43 kDa
Gene Name:	gap junction protein, alpha 1
Database Link:	NP_036699 Entrez Gene 24392 Rat P08050
Synonyms:	Connexin-43; CX43; DFNB38; GJAL; HSS; ODD; ODDD; ODOD; SDTY3



[View online »](#)

Product images:



Western blot of rat cerebellar lysate showing specific immunolabeling of the ~43k connexin43 protein.