

## Product datasheet for **TA309146**

### Slc6a13 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	IHC: 1:100-1:200
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to amino acid residues from the C-terminal region conjugated to KLH
Formulation:	100 µl in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml BSA and 50% glycerol.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	82 kDa
Gene Name:	solute carrier family 6 member 13
Database Link:	<a href="#">NP_598307</a> <a href="#">Entrez Gene 171163 Rat</a> <a href="#">P31646</a>
Synonyms:	GAT-2; GAT2



[View online »](#)