

Product datasheet for **TA309149**

Gfap Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:5,000, IF: 1:1,000
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant and purified bovine GFAP
Formulation:	neat serum
Purification:	Neat Serum
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	50 kDa
Gene Name:	glial fibrillary acidic protein
Database Link:	NP_058705 Entrez Gene 24387 Rat P47819
Synonyms:	FLJ45472



[View online »](#)

Product images:



Western blot of rat cortex lysate showing specific immunolabeling of the ~50k GFAP protein.