

Product datasheet for **TA309204**

Olig1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:3,000, IHC: 1:1000, ICC: 1:1,000, IP: 1:50
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant mouse Olig1
Formulation:	100 µl in 100 mM Glycine, 200 mM Tris, pH 8.0
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	27 kDa
Gene Name:	oligodendrocyte transcription factor 1
Database Link:	NP_058664 Entrez Gene 60394 Rat Entrez Gene 50914 Mouse Q9JKN5
Synonyms:	BHLHB6; BHLHE21; Oligo1



[View online »](#)

Product images:



Western blot of neonatal rat brain lysate showing specific immunolabeling of the ~27k Olig1 protein.