

Product datasheet for **TA309212**

Periostin (POSTN) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000, IHC: 1:100
Reactivity:	Chicken, Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to amino acid residues from the fasciclin domain 1 of mouse periostin
Formulation:	100 ul in PBS
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	93 kDa
Gene Name:	periostin
Database Link:	NP_001129407 Entrez Gene 50706 Mouse Entrez Gene 10631 Human Q15063
Synonyms:	OSF-2; OSF2; PDLPOSTN; PN
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

Product images:



Western blot of rat lung lysate showing specific immunolabeling of the ~ 93 kDa periostin protein triplet.