

Product datasheet for **TA309245**

Synj1 Mouse Monoclonal Antibody [Clone ID: 5H1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	5H1
Applications:	WB
Recommended Dilution:	WB: 1:1000, IHC: 1:100
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein from the C-terminal region of rat synaptojanin 1.
Formulation:	100 µl in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml BSA and 50% glycerol.
Purification:	Protein G purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	145 kDa
Gene Name:	synaptojanin 1
Database Link:	NP_445928 Entrez Gene 85238 Rat Q62910
Synonyms:	INPP5G



[View online »](#)

Product images:



Western blot of rat hippocampal lysate showing specific immunolabeling of the ~145k synaptojanin I protein.