

Product datasheet for **TA309265**

Tubb3 Mouse Monoclonal Antibody [Clone ID: AA10]


Product data:

Product Type:	Primary Antibodies
Clone Name:	AA10
Applications:	WB
Recommended Dilution:	WB: 1:10000, IHC: 1:1000
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Synthetic peptide corresponding to amino acid residues specific to beta III tubulin conjugated to KLH
Formulation:	100 µl in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml BSA and 50% glycerol.
Purification:	Protein G purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55 kDa
Gene Name:	tubulin, beta 3 class III
Database Link:	NP_640347 Entrez Gene 246118 Rat Q4QRB4
Synonyms:	beta-4; MC1R; TUBB4



[View online »](#)

Product images:



Western blot of rat cortex lysate showing specific immunolabeling of the ~55k beta III tubulin protein.