

Product datasheet for **TA309273**

Slc32a1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to amino acid residues from the N-terminal region conjugated to KLH
Formulation:	100 µl in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml BSA and 50% glycerol.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	53 kDa
Gene Name:	solute carrier family 32 member 1
Database Link:	NP_113970 Entrez Gene 83612 Rat O35458
Synonyms:	bA122O1.1; hVIAAT; VGAT; VIAAT



[View online »](#)

Product images:



Western blot of rat hippocampal lysate showing specific immunolabeling of the ~53k VGAT protein.