

Product datasheet for **TA309476**


SHANK3 Mouse Monoclonal Antibody [Clone ID: S69-46]

Product data:

Product Type:	Primary Antibodies
Clone Name:	S69-46
Applications:	WB
Recommended Dilution:	WB, IHC, IP
Reactivity:	Mouse, Human, Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Synthetic peptide amino acids 840-857 of rat Shank3, accession number Q9JLU4
Formulation:	PBS pH7.4, 50% glycerol, 0.09% sodium azide
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	SH3 and multiple ankyrin repeat domains 3
Database Link:	NP_001073889 Entrez Gene 58234 Mouse Entrez Gene 59312 Rat Entrez Gene 85358 Human
Synonyms:	DEL22q13.3; KIAA1650; PROSAP2; PSAP2; SPANK-2
Note:	Detects ~190kDa. No cross-reactivity against Shank1 or Shank2.



[View online »](#)

Product images:

Western blot analysis of Shank3 in rat brain membrane lysates using a 1:1000 dilution of TA309476.