

Product datasheet for **TA309552**

Serotonin N acetyltransferase (AANAT) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	WB, IHC
Reactivity:	Human, Primate
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Human AANAT
Formulation:	Eluted in 0.1N glacial acetic acid containing 0.1% BSA and neutralized using 10x PBS and 0.05% Thimerosal
Purification:	Peptide affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	aralkylamine N-acetyltransferase
Database Link:	NP_001079 Entrez Gene 15 Human Q16613
Synonyms:	DSPS; SNAT
Protein Pathways:	Metabolic pathways, Tryptophan metabolism



[View online »](#)