

## Product datasheet for **TA309570**

### AKAP8 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB, IF, ICC, IHC
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	KLH-conjugated human AKAP95 peptide within amino acids 650 to 692 [Swiss-Prot# O43823].
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	A-kinase anchoring protein 8
Database Link:	<a href="#">NP_005849</a> <a href="#">Entrez Gene 56399 Mouse</a> <a href="#">Entrez Gene 10270 Human</a> <a href="#">O43823</a>
Synonyms:	AKAP-8; AKAP-95; AKAP 95; AKAP95
Protein Families:	Druggable Genome



[View online »](#)