

Product datasheet for **TA309575**


BRCAA1 (ARID4B) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Reacts with residues 1046-1055 [SSKKQKRSHK] of the 136 kDa human, mouse and rat BRCAA1 protein.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	IgG purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	AT-rich interaction domain 4B
Database Link:	NP_057458 Entrez Gene 51742 Human Q4LE39
Synonyms:	BCAA; BRCAA1; RBBP1L1; RBP1L1; SAP180



[View online »](#)

Product images:

Detection of BRCAA1 in MCF-7 whole cell lysate (15 ug) using 1:1,000 dilution of anti-BRCAA1.