

Product datasheet for **TA309576**


Artemin (ARTN) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Amino acids GPGSRARAAGARGC of human Artemin.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	artemin
Database Link:	NP_003967 Entrez Gene 9048 Human Q5T4W7
Synonyms:	artemin; ENOVIN; EVN; NBN; neublastin; neurotrophic factor; neurotrophic factor artemin; OTTHUMP00000009173; OTTHUMP00000009174
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

Product images:

Western Blot: Artemin Antibody - Western Blot analysis of Artemin in human dorsal root ganglia, using anti-Artemin.