

Product datasheet for **TA309578**

ATOX1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ICC/IF, WB
Recommended Dilution:	Immunocytochemistry/ Immunofluorescence: 1:1000, Western Blot: 1:1000
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Full length human ATOX1 protein [Swiss-Prot# O00244]
Formulation:	PBS and 0.1% Sodium Azide
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	antioxidant 1 copper chaperone
Database Link:	NP_004036 Entrez Gene 475 Human O00244
Synonyms:	ATX1; HAH1




[View online »](#)

Product images:



Immunocytochemistry/Immunofluorescence:
ATOX1 Antibody TA309578 - ATOX1 antibody was tested in HeLa cells with Dylight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and Dylight 550 (red).



Western Blot: ATOX1 Antibody TA309578 - HeLa whole cell extract.