

## Product datasheet for **TA309584**

### **NABC1 (BCAS1) Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Reacts with residues 532-550 [DTNSLQNGDKLQKRPEKRQ] of human NABC1.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	IgG purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	breast carcinoma amplified sequence 1
Database Link:	<a href="#">NP_003648</a> <a href="#">Entrez Gene 8537 Human</a> <a href="#">O75363</a>
Synonyms:	AIBC1; NABC1



[View online »](#)

## Product images:



Detection of NABC1 in MCF-7 whole cell lysate (15 ug) using a 1:1,000 dilution of anti-NABC1.