

Product datasheet for **TA309625**

DIABLO Mouse Monoclonal Antibody [Clone ID: 78-1-118]


Product data:

Product Type:	Primary Antibodies
Clone Name:	78-1-118
Applications:	ICC/IF, IP, WB
Recommended Dilution:	Immunocytochemistry/ Immunofluorescence, Western Blot: 1:500, Immunoprecipitation: 1:10 - 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG3/IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant mature human Smac/DIABLO (lacks amino acids 1-55 present in Smac/DIABLO precursor).
Formulation:	Ascitic fluid and 0.1% Sodium Azide
Concentration:	lot specific
Purification:	Ascites
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	diablo IAP-binding mitochondrial protein
Database Link:	NP_063940 Entrez Gene 66593 Mouse Entrez Gene 288753 Rat Entrez Gene 56616 Human Q9NR28
Synonyms:	DFNA64; SMAC
Protein Families:	Transmembrane



[View online »](#)

Product images:



Western Blot: Smac Antibody (78-1-118)
TA309625 - HeLa whole cell extract lysates.