

## Product datasheet for **TA309659**


### **FBXO17 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A synthetic peptide containing residues 70-85 (LQLARDRSAEGRALYA) of the human FBG4 protein. There is 100% sequence identity between the immunogen and mouse and rat sequences.
Formulation:	PBS with 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	F-box protein 17
Database Link:	<a href="#">NP_680474</a> <a href="#">Entrez Gene 50760 Mouse</a> <a href="#">Entrez Gene 292757 Rat</a> <a href="#">Entrez Gene 115290 Human</a> <a href="#">Q96EF6</a>
Synonyms:	FBG4; Fbx17; FBX26; FBXO26
Protein Families:	Druggable Genome



[View online »](#)

**Product images:**

Detection of 30 ug of recombinant human full-length FBG4 using anti-FBG4.