


## Product datasheet for **TA309660**

### Kindlin (FERMT1) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Reacts with residues 643-660 (VHEYIGGYIFLSTRSKDQ) of the 77 kDa human URP1 protein.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	fermitin family member 1
Database Link:	<a href="#">NP_060141</a> <a href="#">Entrez Gene 55612 Human</a> <a href="#">Q9BQL6</a>
Synonyms:	C20orf42; DTGCU2; KIND1; UNC112A; URP1

#### Product images:



Detection of URP1 in 20 ug of human lung adenocarcinoma cell line NCI-H322 lysate using anti-URP1.



[View online »](#)