

## Product datasheet for **TA309678**

### Glutathione Peroxidase 3 (GPX3) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ICC/IF, IHC, Immunoblotting, WB
Recommended Dilution:	Immunohistochemistry-Paraffin: 1:200, Immunoblotting, Immunocytochemistry/ Immunofluorescence: 1:200, Western Blot: 1:500, Immunohistochemistry: 1:200
Reactivity:	Human, Primate, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Full length human GPX3 protein [Swiss-Prot# P22352] expressed in E. coli.
Formulation:	PBS and 0.1% Sodium Azide
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	glutathione peroxidase 3
Database Link:	<a href="#">NP_002075</a> <a href="#">Entrez Gene 64317 Rat</a> <a href="#">Entrez Gene 2878 Human</a> <a href="#">P22352</a>
Synonyms:	GPx-P; GSHPx-3; GSHPx-P
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Arachidonic acid metabolism, Glutathione metabolism




[View online »](#)

Product images:



Immunocytochemistry/Immunofluorescence: Glutathione Peroxidase 3/GPX3 Antibody TA309678 - The GPX3 antibody was tested in HeLa cells at a 1:250 dilution against Dylight 488 (Green). Alpha-tubulin and nuclei were counterstained against Dylight 550 (Red) and DAPI blue).



Western Blot: Glutathione Peroxidase 3/GPX3 Antibody TA309678 - Analysis of GPX3 in human kidney lysate.