

Product datasheet for **TA309702**

ILVBL Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	ICC/IF (Negative), Western Blot: 1:1000
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A synthetic peptide made to an N-terminal portion of the human ILVBL protein (between residues 350-500) [UniProt A1L0T0]
Formulation:	PBS, 30% glycerol and 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ilvB (bacterial acetolactate synthase)-like
Database Link:	NP_006835 Entrez Gene 10994 Human A1L0T0
Synonyms:	209L8; AHAS; ILV2H
Protein Families:	Transmembrane



[View online »](#)

Product images:

250>
150>
100>
75>
50>
37>
25>
20>
15>
10>



Western Blot: ILVBL Antibody TA309702 - Analysis of ILVBL in overexpression lysate.

ILVBL

<250
<150
<100
<75
<50
<37
<25
<20
<15
<10



Western Blot: ILVBL Antibody TA309702 - Analysis of ILVBL in human kidney tissue lysate.