

Product datasheet for **TA309721**

Kallikrein 5 (KLK5) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB, IHC, IP, ELISA
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Residues 280-293 [KFTKWIQETIQANS] of the human KLK-L2 protein.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	kallikrein related peptidase 5
Database Link:	NP_036559 Entrez Gene 25818 Human Q9Y337
Synonyms:	KLK-L2; KLKL2; SCTE
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane

Product images:



Western Blot: KLK L2 Antibody - Western blot of KLK5 (also known as KLK-L2) in BT-474 cells using anti-KLK5 at a 1:1,000 dilution



[View online »](#)