

## Product datasheet for **TA309753**

### MUC4 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB, IF, ICC, IHC
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A genomic peptide made to an internal region of the human MUC4 protein (within residues 1950-2100). [Swiss-Prot Q99102]
Formulation:	PBS, 0.1% BSA, 50% glycerol and 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mucin 4, cell surface associated
Database Link:	<a href="#">NP_004523</a> <a href="#">Entrez Gene 140474 Mouse</a> <a href="#">Entrez Gene 4585 Human</a> <a href="#">Q99102</a>
Synonyms:	ASGP; HSA276359; MUC-4
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)