

Product datasheet for **TA309796**

PDI (PDIA2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ICC/IF, IHC, IP, WB
Recommended Dilution:	Immunohistochemistry: 1:100-1:1000, Immunocytochemistry/ Immunofluorescence: 1:5000-1:20000, Immunoprecipitation: 1-3 ug, Western Blot: 1:10000-1:20000, Immunohistochemistry-Paraffin: 1:100-1:500, Immunohistochemistry-Frozen: 1:100-1:1000
Reactivity:	Bovine, Human, Mouse, Primate, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Purified full-length liver protein.
Concentration:	lot specific
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	protein disulfide isomerase family A member 2
Database Link:	NP_006840 Entrez Gene 69191 Mouse Entrez Gene 287164 Rat Entrez Gene 64714 Human Q13087
Synonyms:	PDA2; PDI; PDIP; PDIR
Protein Families:	Druggable Genome



[View online »](#)