

## Product datasheet for **TA309809**

### **PMP22 Mouse Monoclonal Antibody [Clone ID: Hu1]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	Hu1
Applications:	IHC
Recommended Dilution:	IHC
Reactivity:	Human, Primate
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Human PMP22 cDNA boosted with a 13-amino acid peptide of the 2nd extracellular domain (a.a. 120-133) conjugated to murine IgG. [Swiss-Prot# Q01453]
Formulation:	Tris-glycine, 150mM NaCl and 0.05% Sodium Azide
Concentration:	lot specific
Purification:	IgG purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	peripheral myelin protein 22
Database Link:	<a href="#">NP_000295</a> <a href="#">Entrez Gene 5376 Human Q01453</a>
Synonyms:	CMT1A; CMT1E; DSS; GAS-3; GAS3; HMSNIA; HNPP; Sp110
Protein Families:	Transmembrane



[View online »](#)