

Product datasheet for **TA309813**

DNA Polymerase iota (POLI) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A 15-mer peptide in the C-terminal region of human DNA Polymerase Iota, conjugated to KLH.
Formulation:	PBS and 0.02% Sodium Azide
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	polymerase (DNA) iota
Database Link:	NP_009126 Entrez Gene 26447 Mouse Entrez Gene 11201 Human Q9UNA4
Synonyms:	RAD30B; RAD30B
Protein Families:	Druggable Genome



[View online »](#)