

Product datasheet for **TA309821**


PTCH1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB, IHC, ICC, ELISA
Reactivity:	Chicken, Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A synthetic peptide made to residues 269-279 of the human and mouse PTCH protein.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	IgG purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	patched 1
Database Link:	NP_000255 Entrez Gene 19206 Mouse Entrez Gene 89830 Rat Entrez Gene 5727 Human Q13635
Synonyms:	BCNS; HPE7; NBCCS; PTC; PTC1; PTCH; PTCH11
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Basal cell carcinoma, Hedgehog signaling pathway, Pathways in cancer



[View online »](#)

Product images:

Western Blot: Patched 1 Antibody -Analysis of Human PTCH1 in Fetal brain Lysate with anti huPTCH1 at a 1:500 dilution.