

Product datasheet for **TA309822**


PTCH2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB, IHC, ICC, ELISA
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	residues 226-235 [GPFASLEGFR] of the human PTCH2 protein
Formulation:	0.1 mg in 0.1 ml (1mg/1ml) in PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	patched 2
Database Link:	NP_001159764 Entrez Gene 19207 Mouse Entrez Gene 8643 Human Q9Y6C5
Synonyms:	PTC2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Basal cell carcinoma, Hedgehog signaling pathway, Pathways in cancer



[View online »](#)

Product images:

Western Blot: Patched 2 Antibody -Analysis of Human PTCH2 in Fetal brain Lysate with Anti huPTCH2 at a 1:500 dilution.