

Product datasheet for **TA309834**

Calcipressin 1 (RCAN1) Rabbit Polyclonal Antibody


Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | ICC/IF, IHC, WB |
| Recommended Dilution: | Immunohistochemistry: 1:1000, Western Blot: 2 ug/ml, Immunohistochemistry-Paraffin: 1:1000, Immunocytochemistry/ Immunofluorescence |
| Reactivity: | Human, Mouse, Rat |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic peptide made to an internal portion of the human protein (within residues 210-252). [Swiss-Prot# P53805] |
| Formulation: | PBS with 30% glycerol and 0.1% Sodium Azide |
| Concentration: | lot specific |
| Purification: | Affinity purified |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Gene Name: | regulator of calcineurin 1 |
| Database Link: | NP_981963 Entrez Gene 54720 Mouse Entrez Gene 266766 Rat Entrez Gene 1827 Human P53805 |
| Synonyms: | ADAPT78; CSP1; DSC1; DSCR1; MCIP1; RCN1 |
| Protein Families: | Transcription Factors |




[View online »](#)

Product images:



Immunohistochemistry: RCAN1 Antibody TA309834 - Analysis of RCAN1 in human fetal brain



Western Blot: RCAN1 Antibody TA309834 - Analysis of RCAN in human fetal brain