


Product datasheet for TA309837

RFXANK Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human, Primate
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Peptide made to a portion of the RFX-B protein.
Formulation:	PBS. and 0.02% Sodium Azide
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	regulatory factor X associated ankyrin containing protein
Database Link:	NP_003712 Entrez Gene 8625 Human O14593
Synonyms:	ANKRA1; BLS; F14150_1; RFX-B
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Antigen processing and presentation, Primary immunodeficiency

Product images:



Western Blot: RFXANK Antibody - Western Blot analysis of RFXANK on Raji lymphoblastoid cell line at 1:500 dilution.



[View online >](#)