

Product datasheet for **TA309841**


RPS27A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF
Recommended Dilution:	WB, IF, IHC, ELISA
Reactivity:	Bovine, C. elegans, Chicken, Drosophila, Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Gluteraldehyde cross-linked ubiquitin.
Purification:	Whole antisera
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ribosomal protein S27a
Database Link:	NP_002945 Entrez Gene 78294 Mouse Entrez Gene 6233 Human P62979
Synonyms:	CEP80; HEL112; S27A; UBA80; UBC; UBCEP1; UBCEP80
Protein Families:	Druggable Genome
Protein Pathways:	Ribosome



[View online »](#)

Product images:

Immunofluorescent: Ubiquitin Antibody -
Immunofluorescent staining of aggregates of
mouse L-cells treated with lactacystine, using
anti-Ubiquitin.