

## **Product datasheet for TA309842**

## OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **RUNX2 Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

**Applications:** FC, ICC/IF, IHC

Recommended Dilution: Immunohistochemistry-Paraffin: 1:200, Immunohistochemistry: 1:200,

Immunocytochemistry/ Immunofluorescence: 1:50-1:100, Flow Cytometry: 2-5 ug/ml

Reactivity: Human, Mouse

**Host:** Rabbit

Clonality: Polyclonal

**Immunogen:** A synthetic peptide made to an internal region of the human RUNX2 protein (within residues

225-300). [Swiss-Prot Q13950]

**Formulation:** PBS, 30% glycerol and 0.05% Sodium Azide

Concentration: lot specific

Purification: Affinity purified
Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** runt related transcription factor 2

Database Link: NP 004339

Entrez Gene 12393 MouseEntrez Gene 860 Human


Q13950

Synonyms: AML3; CBFA1; CCD; CCD1; OSF-2; OSF2; PEA2aA; PEBP2A1; PEBP2A2; PEBP2aA; PEBP2aA1


**Protein Families:** Druggable Genome, Transcription Factors



## **Product images:**




Immunocytochemistry/Immunofluorescence: RUNX2/CBFA1 Antibody TA309842 - RUNX2 antibody was tested in HeLa cells with FITC (green). Nuclei and alpha-tubulin were counterstained with Dapi (blue) and Dylight 550 (red).



Flow Cytometry: RUNX2/CBFA1 Antibody TA309842 - An intracellular stain was performed on HeLa cells with RUNX2/CBFA1 Antibody TA309842AF647 (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 5 ug/mL for 30 minutes at room temperature. Both antibodies were conjugated to Alexa Fluor 647.





Immunohistochemistry: RUNX2/CBFA1 Antibody TA309842 - Staining of RUNX2 in mouse prostate using DAB with hematoxylin counterstain.