

Product datasheet for **TA309854**

SF3B4 Mouse Monoclonal Antibody [Clone ID: 3A1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3A1
Applications:	WB
Recommended Dilution:	WB, IF
Reactivity:	Bovine, Human, Mouse, Porcine, Rat
Host:	Mouse
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Immunogen:	Full length recombinant human SF3B4 expressed in and purified from E. coli.
Formulation:	PBS and 10mM Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	splicing factor 3b subunit 4
Database Link:	NP_005841 Entrez Gene 107701 MouseEntrez Gene 295270 RatEntrez Gene 10262 Human Q15427
Synonyms:	AFD1; Hsh49; SAP49; SF3b49
Protein Pathways:	Spliceosome



[View online »](#)

Product images:

Western Blot: SF3B4 Antibody (3A1) - Western blot of HeLa cell crude extract stained with anti-SF3B4. The SF3B4 (3A1) antibody runs with an apparent SDS-PAGE molecular weight 49kDa.