

## Product datasheet for **TA309857**


### SHANK1 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IF, IHC
Recommended Dilution:	WB, IHC, IF
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	residues [SGPIYPGLFDIRSS] of the C-terminus of the Shank1a protein.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	SH3 and multiple ankyrin repeat domains 1
Database Link:	<a href="#">NP_057232</a> <a href="#">Entrez Gene 78957 Rat</a> <a href="#">Entrez Gene 243961 Mouse</a> <a href="#">Entrez Gene 50944 Human</a> <a href="#">Q9Y566</a>
Synonyms:	SPANK-1; SSTRIP; synamon



[View online »](#)

**Product images:**

Immunohistochemistry: SHANK1a Antibody -  
Immunohistochemical staining of mouse retina  
using anti-SHANK1a.



Mouse C57B6 staining with anti-SHANK1a.