

Product datasheet for **TA309858**

SHANK1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide: ERLIREYPQSFKEGV, corresponding to C terminal amino acids 135-149 of Human SHANK1.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	SH3 and multiple ankyrin repeat domains 1
Database Link:	NP_057232 Entrez Gene 78957 Rat Entrez Gene 243961 Mouse Entrez Gene 50944 Human Q9Y566
Synonyms:	SPANK-1; SSTRIP; synamon

Product images:



Western Blot: SHANK1 Antibody - Analysis of SHANK1 in Rat Brain lysate at a 1:500 dilution.



[View online »](#)