

Product datasheet for **TA309860**

SIRT7 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IP, WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide made to an internal portion of the human protein (within residues 100-200). [Swiss-Prot# Q9NRC8]
Formulation:	Tris-citrate/phosphate, pH 7-8 and 0.1% Sodium Azide
Concentration:	lot specific
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	sirtuin 7
Database Link:	NP_057622 Entrez Gene 51547 Human Q9NRC8
Synonyms:	SIR2L7
Protein Families:	Druggable Genome, Transcription Factors




[View online »](#)

Product images:



Detection of SIRT7 in human liver lysates using TA309860.



SIRT7 Immunoprecipitated from human liver lysates with A). [TA309859] B). TA309860