

Product datasheet for **TA309861**

VMAT2 (SLC18A2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB, IHC
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A genomic peptide made to an internal region of the human VMAT2 protein (within residues 30-200). [Swiss-Prot Q05940]
Formulation:	PBS, 30% glycerol and 0.05% Sodium Azide
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	solute carrier family 18 member A2
Database Link:	NP_003045 Entrez Gene 25549 Rat Entrez Gene 214084 Mouse Entrez Gene 6571 Human Q05940
Synonyms:	SVAT; SVMT; VAT2; VMAT2
Protein Families:	Transmembrane
Protein Pathways:	Parkinson's disease



[View online »](#)