

## Product datasheet for **TA309932**

### WWOX Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	FC, IHC
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A synthetic peptide corresponding to residues 89-107 [RLAFTVDDNPTKPTTRQRY] of the human protein.
Formulation:	PBS and 0.02% Sodium Azide
Concentration:	lot specific
Purification:	IgG purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	WW domain containing oxidoreductase
Database Link:	<a href="#">NP_057457</a> <a href="#">Entrez Gene 80707 Mouse</a> <a href="#">Entrez Gene 51741 Human</a> <a href="#">Q9NZC7</a>
Synonyms:	D16S432E; FOR; FRA16D; HHCMA56; PRO0128; SCAR12; SDR41C1; WOX1
Protein Families:	Druggable Genome



[View online »](#)