

Product datasheet for **TA309937**


Huntingtin (HTT) Mouse Monoclonal Antibody [Clone ID: H.F5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	H.F5
Applications:	WB
Recommended Dilution:	WB, IHC
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	N-terminal peptides 1-225 aa.
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	huntingtin
Database Link:	NP_002102 Entrez Gene 15194 Mouse Entrez Gene 3064 Human P42858
Synonyms:	HD; IT15
Protein Families:	Druggable Genome
Protein Pathways:	Huntington's disease



[View online »](#)

Product images:

Western blot analysis of Huntingtin expression in HEK293 cells.