

Product datasheet for **TA310046**


Ptgs1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Murine COX 1.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	prostaglandin-endoperoxide synthase 1
Database Link:	NP_032995 Entrez Gene 24693 Rat Entrez Gene 19224 Mouse P22437
Synonyms:	COX-1; COX1; COX3; Cyclooxygenase-1; PCOX1; PGG/HS; PGHS-1; PGHS1; PHS1; PTGHS



[View online »](#)

Product images:

Western blot analysis of Jurkat cells (Lanes 1 and 2: 62ug cell lysate) with anti-COX1 antibody. (Predicted band size: 70kDa; observed band size: 50kDa)