

## Product datasheet for **TA310060**

### Kallikrein 14 (KLK14) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:100
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Purification:	Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	kallikrein related peptidase 14
Database Link:	<a href="#">NP_071329</a> <a href="#">Entrez Gene 43847 Human</a> <a href="#">Q9P0G3</a>
Synonyms:	KLK-L6
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)

**Product images:**

Western blot analysis of fetal liver lysates with anti-Kallikrein 14 antibody (1:100).