

Product datasheet for **TA310066**


MTH1 (NUDT1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide conjugated to KLH
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Purification:	Serum
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	nudix hydrolase 1
Database Link:	NP_002443 Entrez Gene 4521 Human P36639
Synonyms:	MTH1
Protein Families:	Stem cell - Pluripotency



[View online »](#)

Product images:

Western blot analysis of HEK293 and HeLa cell lysates with anti-MTH1 antibody (1:500). This antibody identified a band of ~18kDa.