

# Product datasheet for TA310076

## PGC1 alpha (PPARGC1A) Rabbit Polyclonal Antibody

### **Product data:**

#### OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Human PGC-1.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	PPARG coactivator 1 alpha
Database Link:	<u>NP_037393</u> <u>Entrez Gene 19017 MouseEntrez Gene 83516 RatEntrez Gene 10891 Human</u> <u>Q9UBK2</u>
Synonyms:	LEM6; PGC-1(alpha); PGC-1alpha; PGC-1v; PGC1; PGC1A; PPARGC1
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Adipocytokine signaling pathway, Huntington's disease, Insulin signaling pathway



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

#### **Product images:**



Western blot analysis of Thp-1 cells (Lane 1: 62ug cell lysate) and ECL cells (Lane 2: 70ug cell lysate) with anti-PGC1 antibody diluated at 1:1000. (Predicted band size: 89kDa; observed band size: 55kDa)




Figure from citation: Western blot analysis of PGC-1a (also known as PPARGC1A) protein level by using anti-PPARGC1A antibody in control and PGC-1a-overexpressing H9c2 cardiomyocytes. <u>View Citation</u>

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US