

Product datasheet for **TA310083**


Vesicle docking protein p115 (USO1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:250
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide conjugated to KLH
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Purification:	Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	USO1 vesicle transport factor
Database Link:	NP_003706 Entrez Gene 8615 Human O60763
Synonyms:	P115; TAP; VDP



[View online »](#)

Product images:

Western blot analysis of fetal kidney lysates with anti-VDP antibody (1:250). This antibody identified a band of ~110kDa in human fetal kidney lysates.