

Product datasheet for **TA310087**

ABHD6 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Human ABHD6.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	abhydrolase domain containing 6
Database Link:	NP_065727 Entrez Gene 66082 Mouse Entrez Gene 305795 Rat Entrez Gene 57406 Human Q9BV23
Synonyms:	abhydrolase domain containing 6; lipase protein
Protein Families:	Transmembrane



[View online »](#)

Product images:

Western blot analysis of THP1 cells (Lane 1: 60ug lysate) and ECL cells (Lane 2: 72ug lysate) with anti-ABHD6 antibody diluted at 1:2000. (Predicted and observed size: 38kDa)