

Product datasheet for **TA310092**


ADAMTS2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide conjugated to KLH
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Purification:	Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ADAM metallopeptidase with thrombospondin type 1 motif 2
Database Link:	NP_055059 Entrez Gene 9509 Human O95450
Synonyms:	ADAM-TS2; ADAMTS-2; ADAMTS-3; NPI; PC I-NP; PCI-NP; PCINP; PCPNI; PNPI
Protein Families:	Druggable Genome



[View online »](#)

Product images:

Western blot analysis of human fetal brain lysates with anti-ADAMTS2 antibody (1:1000). This antibody identified zymogen form of ADAMTS2 at 130-140kDa, activated forms at 94-105kDa, and breakdown products at 50kDa and 34kDa in human fetal brain lysates.