

## Product datasheet for **TA310093**


### ADAMTS3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:2000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide conjugated to KLH
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Purification:	Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ADAM metallopeptidase with thrombospondin type 1 motif 3
Database Link:	<a href="#">NP_055058</a> <a href="#">Entrez Gene 9508 Human</a> <a href="#">O15072</a>
Synonyms:	ADAMTS-4
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)

**Product images:**

Western blot analysis of human fetal heart (Lane 1) and brain (Lane 2) lysates with anti-ADAMTS3 antibody (1:2000). This antibody identified zymogen form of ADAMTS3 at 140kDa, and activated forms at 94-105kDa.