

Product datasheet for **TA310094**


ADAMTS5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide conjugated to KLH
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Purification:	Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ADAM metallopeptidase with thrombospondin type 1 motif 5
Database Link:	NP_008969 Entrez Gene 11096 Human Q9UNA0
Synonyms:	ADAM-TS 5; ADAM-TS5; ADAM-TS 11; ADAMTS-5; ADAMTS-11; ADAMTS11; ADMP-2
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)

Product images:

Western blot analysis of human fetal brain lysates with anti-ADAMTS5 antibody (1:500). This antibody identified a band of 72kDa and a breakdown product at 50kDa in human fetal brain tissue lysates.