

Product datasheet for **TA310101**

AGMAT Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:200~1000, IHC: 1:100~500, ELISA: 1:5000~10,000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Rabbit polyclonal AGMAT antibody was raised against a peptide from the intermediate residues of human AGMAT protein.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	agmatinase
Database Link:	NP_079034 Entrez Gene 79814 Human Q9BSE5
Synonyms:	FLJ23384
Protein Pathways:	Arginine and proline metabolism, Metabolic pathways



[View online »](#)

Product images:



Western blot analysis of fetal kidney lysates with anti-AGMAT antibody (1:500).