

Product datasheet for **TA310112**


Dicer (DICER1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 1:100~1:500
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide conjugated to KLH
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	dicer 1, ribonuclease III
Database Link:	NP_803187 Entrez Gene 192119 Mouse Entrez Gene 23405 Human Q9UPY3
Synonyms:	DCR1; Dicer; Dicer1e; HERNA; K12H4.8-LIKE; MNG1; RMSE2
Protein Families:	Druggable Genome



[View online »](#)

Product images:

Immunohistochemical staining of formalin-fixed paraffin-embedded brain showing cytoplasmic staining with anti-Dicer1 antibody diluted at 1:100.