

Product datasheet for **TA310115**

HCLS1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:200~1000, IHC: 1:50~200, ELISA: 1:40,000~80,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide conjugated to KLH
Formulation:	PBS, pH 7.4 with 0.05% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	hematopoietic cell-specific Lyn substrate 1
Database Link:	NP_005326 Entrez Gene 15163 Mouse Entrez Gene 288077 Rat Entrez Gene 3059 Human P14317
Synonyms:	CTTNL; HS1; lckBP1; p75
Protein Families:	Transcription Factors
Protein Pathways:	Pathogenic Escherichia coli infection, Tight junction



[View online »](#)

Product images:



Western blot analysis of Raji cell lysates with anti-HCLS1 antibody (1:500).